

## Side-by-side Door Style Blood Bank Refrigerator

### Advanced Cooling System

Ensure temperature fast recover and stability by intelligent controlled fan forced air cooling system, air cooled chiller finned evaporator with unique circulation system design to ensure frost-free cooler environment, double-layer vacuum electric heating glass door, ensure the door will not get condensation under high temperature and high humidity environment.

### Precised Microcomputer temperature control

Microprocessor-based temperature controller, digital temperature display, accuracy is  $0.1^{\circ}\text{C}$ . Excellent refrigeration layout, inside temperature is  $4\pm 1^{\circ}\text{C}$ .

Standard thermal printer can track temperature inside all the way.

### Operating Security

For safety and reliability, alarm module keep monitoring surroundings and system status

Two Alarm style: buzzing, flashing

ALarms: High/Low temperature alarms, power failure alarm for 72 hours, sensor failure alarm, door ajar alarm, condenser abnormal alarm, low battery alarm.

### Humanized Design

Equipd baskets for blood bags.

Safety lock design prevents any abnormal door status.

Test hole optional

### Adjustable shelves design.

Lights inside make everything clear.

Foam glass door design, independent inner doors, self-closed doors ensure perfect insulation performance.

Patented temperature box design, make the display temperature closer to the storage blood Sus inner tank

Wide voltage range from 187V~242V.

### Side-by-Side Door Style



## Henan Forever Medical CO., LTD

### Technical Parameters :

Model	BBR660	BBR1000
Cabinet Type	Vertical	
Refrigerant	CFC-Free	
Cooling Performance(Centigrade)	4	
Temperature Range(Centigrade)	4+-1	
Capacity	660L	1000L
Power	364W	513W
Noise (db)	51.2	53.6
Climate	N	
N.W.(kg)	165/195	256/286
Capacity/Baskets/450L Blood Bag/piece	658/36/360	1000/12/576
Internal size(W*D*H/mm)	1100*454*1325	1100*684*1325
External size(W*D*H/mm)	1220*648*1885	1220*878*1885
Package Size(W*D*H/mm)	1300*730*2080	1300*960*2080
Alarms	High and low temperature alarms,sensor error alarm, power failure alarm, door ajar alarm,low battery alarm	
Accessories	6 inner doors,testing hole,out door latch,light,thermal printer	
Certificate	ISO	